

Brain Nuclear Extracts

(2/97, AG; modified from David Frank)

Lysis Procedure

1. Place small piece of frozen ctx (approx one hemisphere of E18 Ctx) in 1 ml PBS lysis buffer.
2. Homogenize with dounce homogenizer (5 strokes with A, 5 strokes with B).
3. Transfer homogenate to eppendorff tube. Spin at 4 C for 5 mins.
4. Discard supernatant, resuspend pellet (nuclei) in 1 ml Extraction Lysis Buffer.
5. Sonicate nuclei with two 15 sec bursts.
6. Spin at 4 C for 15 mins.
7. Recover supernatant and transfer to new tube.
8. Quantify protein concentration by Bradford Assay. Label tubes and freeze at -70 C.
9. To run lysate in gel, suspend required amount in 2X SDS sample buffer (50 μ l) with b-mercapto-ethanol, boil, cool on ice, and load.

Note: To make whole brain lysate, go directly to step 4 and start by placing tissue in 1 ml extraction buffer.

PBS Lysis Buffer (make fresh with cold PBS):

PBS 5 mls

200mM PMSF 25 μ l

2mg/ml Aprotinin 25 μ l

Extraction Buffer (for 100mls; store at 4 C):

ddH₂O 83 mls

1M Tris (pH7.8) 2 mls (20 mM)

5M NaCl 2.5 mls (125 mM)

1M MgCl₂ 0.5 mls (5 mM)

0.5 M EDTA 40 μ l (0.2 mM)

Glycerol 12 mls (12%)

NP40 100 μ l (0.1%)

Just before using prepare **Extraction Lysis Buffer** as follows:

Extraction buffer 5 mls

Aprotinin 25 μ l

PMSF (200 mM) 25 μ l

200 mM DTT 250 μ l

Leupeptin (10 mg/ml) 5 μ l