

Slice overlay assay

[as used in Polleux *et al.* (1998) *Science* 282:1904.]

Cortical slices

Embryonic day 18 (E18) to postnatal day 3 (PND3) rat brains were rapidly isolated in ice-cold HBSS and then embedded in 2.5 percent low-melting point agarose (diluted in HBSS) placed on ice to accelerate solidification. 250 microns thick coronal sections were performed using a vibratome, filled with HBSS.

Sections are collected using a dropper (large opening flamed-tip Pasteur pipette) and plated onto a membrane insert of a 6-wells plate (Becton Dickinson-1 micron pore size). Put 1.8 mls of slice culture medium beneath the insert about 30 minutes before plating and leave the plate in a cell culture incubator (37°C-5 percent CO₂).

Dissociated cells

After enzymatic dissociation of E18 rat cortex (see [dissociation protocol](#)), cells were labeled with the carbocyanine fluorescent dye DiI (Molecular Probes; diluted 10 mg/ml in 100 percent ethanol). Briefly, 10 microliters of the DiI stock solution (ultrasonicate right before using it) is added to 1 ml of cell suspension (5x10⁶ cells/ml) and incubated for 15 minutes @ 37°C. Then the cell suspension is centrifuged gently (for 15 ml conical tubes, 1000 rpm-5 minutes @ room temperature) and the supernatant is then discarded. The pellet is gently resuspended in 5 mls of fresh, prewarmed cell culture medium. This washing step is repeated 3 more times to get rid of DiI microscopic crystals.

After the final wash, cells are resuspended @ a final concentration of 5x10⁵ cells/ml and plated onto cortical slices prepared about an hour before the cells. The plating part can be tricky ! Delicately pipet about 400 microliters of the DiI labeled cells solution onto and around the slice. Don't drop the solution ! Shake the plate in X and Y axis gently three or four times and then put the plate in a cell culture incubator. The axon outgrowth can be scored typically after 3-5 hours using live cell imaging.

Solutions

Hank's Balanced Salt Solution (HBSS)

10X HBSS (Gibco #310-4180) 50 ml

1 M Hepes (pH 7.4) 1.25 ml (=2.5mM)

1 M Glucose 15 ml (=6.5 mg/ml=35 mM)

100 mM CaCl₂ 5 ml (=1mM)

100 mM MgSO₄ 5 ml (=1mM)

1 M NaHCO₃ 2 ml (=4mM)

Add ddH₂O to a total vol. of 500 ml. Filter sterile.

Slice culture medium for E18 to P7 rat brain.

For 50 mls:

34.5 mls Basal Medium Eagle (without L-glutamine)

12.5 mls HBSS (see above)

20 mM Glucose

1 mM L-glutamine

1mM Penicilin-Streptomycin

Filter sterile then add 5 percent Normal Horse Serum (heat inactivated).

For postnatal slice culture add kynurenic acid/Mg²⁺ (broad NMDA antagonist to prevent glutamate induced excitotoxicity) to the medium.

Cell culture medium for rat E18 culture.

See [solutions for cell culture](#).